



A



B



C



D




E



F

EYE INSTRUMENTS

- | | | |
|---|---------|--|
| A | 420-802 | Eye Dressing Forceps, Straight, Extra Fine, Serrated, 4" |
| B | 420-810 | Eye Dressing Forceps, Half Curved, Extra Fine, Serrated, 4" |
| C | 420-828 | Eye Dressing Forceps, Fully Curved, Extra Fine, Serrated, 4" |
| D | 421-800 | Iris Tissue Forceps, Straight, Extra Fine, 1 x 2 Teeth, 4" |
| E | 421-818 | Iris Tissue Forceps, Half Curved, Extra Fine 1 x 2 Teeth, 4" |
| F | 421-826 | Iris Tissue Forceps, Fully Curved, Extra Fine, 1 x 2 Teeth, 4" |



Eye Instruments

- 740-050 Plain Wire Eye Speculum, 2-1/2" (6.4 cm), blades 15 mm wide
- 740-506 DESMARRES Lid Retractor, 5-1/2" (14 cm), size 1, 11 mm
- 740-514 DESMARRES Lid Retractor, 5-1/2" (14 cm), size 2, 13 mm
- 740-522 DESMARRES Lid Retractor, 5-1/2" (14 cm), size 3, 15 mm
- 740-524 DESMARRES Lid Retractor, 5-1/2" (14 cm), size 4, 17 mm
- 741-009 GRAEFE Eye Dressing Forceps, 2-3/4" (7 cm), straight, serrated
- 741-017 GRAEFE Eye Dressing Forceps, 2-3/4" (7 cm), curved, serrated
- 741-157 GRAEFE Iris Forceps, 2-3/4" (7 cm), straight, 1 X 2 teeth
- 741-165 GRAEFE Iris Forceps, 2-3/4" (7 cm), curved, 1 X 2 teeth



EYE INSTRUMENTS

- A** 421-830 Foerster Eye Dressing Forceps, 3-3/4", Octagonal Grip, Straight, Serrated Tips, 0.5mm Wide
421-832 Foerster Eye Dressing Forceps, 3-3/4", Octagonal Grip, Straight, Smooth Tips, 0.5mm Wide
- B** 421-835 Foerster Eye Dressing Forceps, 3-3/4", Octagonal Grip, Curved, Serrated Tips, 0.5mm Wide
421-838 Foerster Eye Dressing Forceps, 3-3/4", Octagonal Grip, Curved, Smooth Tips, 0.5mm Wide
- C** 421-840 Foerster Iris Tissue Forceps, 3-3/4", Octagonal Grip, Straight, 1 x 2 Teeth, 0.5mm Wide
- D** 421-845 Foerster Iris Tissue Forceps, 3-3/4", Octagonal Grip, Curved, 1 x 2 Teeth, 0.5mm Wide



CASTROVIEJO

A



B



C



WILLIAMS 2 3/4" (7 cm)
3 1/4" (8.2 cm)

D



GRAEFE 3 1/2" (8.9 cm)

E



LANCASTER

F



G

EYE INSTRUMENTS SPECULUMS

- A 741-220 Castroviejo Eye Speculum, 3.75" Small Blades 13x5mm O.D
- 741-222 Castroviejo Eye Speculum, 3.75" Medium Blades 15x5mm O.D.
- 741-224 Castroviejo Eye Speculum, 4" Large Blades 18x5mm O.D.

- B 741-250 Murdock Eye Speculum, 2.25" Blades 15mm Wide


- C 741-300 Barraquer Eye Speculum, 1.75" Small Blades 10mm Wide
- 741-310 Barraquer Eye Speculum, 1.75" Large Blades 15mm Wide

- D 741-320 Williams Eye Speculum 2 3/4", Small Blades 11x5mm
- 741-325 Williams Eye Speculum 3 3/4" Large Blades 15x5mm

- E 741-350 Von Graefe Eye Speculum, 3.5"


- F 741-355 Lancaster Eye Speculum, 3.25" Blades 15x5mm

- G 741-360 Guyton-Park Eye Speculum 3.5" Solid Blades 14mm Wide with Suture Posts
- 741-365 Guyton-Park Eye Speculum, 3.5" Fenestrated Blades 14mm Wide with Suture Posts



EYE INSTRUMENTS SPECULUMS

- | | | |
|----------|----------------|--|
| A | 741-370 | Barraquer Eye Speculum, 1 ½" Long, 10mm, Solid Blade, Child |
| | 741-372 | Barraquer Eye Speculum, 1 ½" Long, 15mm, Solid Blade, Adult |
| B | 741-375 | Smith Eye Speculum, 5cm |
| C | 741-380 | Clark Eye Speculum, 8cm |
| D | 741-385 | Knapp Eye Speculum, Non-Magnetic |
| E | 741-390 | Mellinger Eye Speculum, 5cm |
| | 741-395 | Mellinger Eye Speculum, 7cm |
| F | 741-400 | Jaeger Lid Plate, 11cm |




EYE INSTRUMENTS

- A 741-504 Von Graefe Fixation Forceps, 4-1/4"
B 741-512 Von Graefe Fixation Forceps, with Catch, 4-1/4"
C 741-524 Bishop-Harmon Dressing Forceps, Str., Del., 0.5mm Wide, 3-3/8"
D 741-528 Bishop-Harmon Dressing Forceps, Straight, Serr., 0.8mm Wide, 3-3/8"
E 741-532 Bishop-Harmon Tissue Forceps, Str., 1 x 2 Teeth, 0.5mm Wide, 3-3/8"
F 741-536 Bishop-Harmon Tissue Forceps, Str., 1 x 2 Teeth, 0.8mm Wide, 3-3/8"



741-603
LESTER Fixation
delicate, 1 X 2 teeth



741-605
LESTER Fixation
delicate, 2x3 teeth



741-801



742-007



742-502



742-551



742-759

EYE INSTRUMENTS

- 741-603 LESTER Fixation Forceps, 3-3/4" (9.5 cm), delicate, 1 X 2 teeth, 0.6 mm wide
- 741-605 LESTER Fixation Forceps, 3-3/4" (9.5 cm), delicate, 2x3 teeth, 1 mm wide
- 741-801 Eye Knife test drum
- 742-007 DESMARRES Chalazion Forceps, 3-1/2" (8.9 cm), large size, inside ring 17 X 28 mm
- 742-502 BARRAQUER Blade Breaker and Holder, 4-1/2" (11.4 cm), with ratchet lock
- 742-551 CASTROVIEJO Blade Breaker and Holder, 5-1/4" (13.3 cm), with push-pull lock
- 742-759 WELLS Enucleation Scoop, 6" (15.2 cm), cup 21 mm



742-908
742-916
742-924



00



743-005
743-013
743-021
743-023
743-025
743-028



743-252



743-302

EYE INSTRUMENTS

- | | |
|---------|--|
| 742-908 | Von Graefe Strabismus Hook, Small, 5 1/2" |
| 742-916 | Von Graefe Strabismus Hook, Medium, 5 1/2" |
| 742-924 | Von Graefe Strabismus Hook, Large, 5 1/2" |
| 743-005 | Meyhoeffler Chalazion Curette, 5", Size 00, 1mm |
| 743-013 | Meyhoeffler Chalazion Curette, 5", Size 0, 1.5mm |
| 743-021 | Meyhoeffler Chalazion Curette, 5", Size 1, 1.7mm |
| 743-023 | Meyhoeffler Chalazion Curette, 5", Size 2, 2mm |
| 743-025 | Meyhoeffler Chalazion Curette, 5", Size 3, 2.5mm |
| 743-028 | Meyhoeffler Chalazion Curette, 5", Size 4, 3.5mm |
| 743-252 | Risley (LaForce) Foreign Body Spud, Angular, 6" |
| 743-302 | Dix Needle and Spud, Reversible Handle, 6" |



A



B



C



D



E




F



G


EYE INSTRUMENTS

- | | | |
|---|---------|--|
| A | 742-928 | Green Strabismus Hook, 5" 9mm Deep |
| B | 742-930 | O'Connor Tenotomy Hook, 5", Fine serrations inside |
| C | 742-934 | Jameson Strabismus Hook, 5", 10mm Deep |
| D | 742-936 | New Orleans Lens Loop, 5 1/4", Smooth Loop, 4.5 x 10mm Very Delicate |
| E | 742-938 | Kirby Muscle Hook and Expressor 6", Double End |
| F | 742-940 | Wilder Lens Loop, 5 1/4", Serrated Loop, 4x15mm |
| G | 742-942 | Wilder Lens Loop, 5 1/4", Smooth Loop, 4x15mm |



EYE INSTRUMENTS

- | | | |
|---|---------|--|
| A | 743-450 | Bowman Probe, 000-0000, Double Ended, Sterling Silver, 5" |
| | 743-468 | Bowman Probe, 00-0, Double Ended, Sterling Silver, 5" |
| | 743-476 | Bowman Probe, 1-2, Double Ended, Sterling Silver, 5" |
| | 743-484 | Bowman Probe, 3-4, Double Ended, Sterling Silver, 5" |
| | 743-492 | Bowman Probe, 5-6, Double Ended, Sterling Silver, 5" |
| | 743-500 | Bowman Probe, 7-8, Double Ended, Sterling Silver, 5" |
| B | 743-505 | Thorpe Caliper 4 1/2" (11.4cm), Graduated in inches and mm, Permits measurements in deep areas |
| C | 743-807 | Wilder Lacrimal Dilator, 4", Fine Taper |
| | 743-815 | Wilder Lacrimal Dilator, 4", Medium Taper |
| | 743-816 | Wilder Lacrimal Dilator, 4", Heavy Taper |
| D | 743-820 | Jameson Caliper 3 3/4" (9.5cm), Graduated in inches and mm, with thin tips, Chrome |



EYE INSTRUMENTS

- A** 103-634 Iris Scissors, 4-1/2" (11.4 cm), angled on side
- B** 744-056 Eye Scissors, Large Rings, Sharp Tips, Straight, 4", Blades (104 mm)
- 744-058 Eye Scissors, Large Rings, Sharp Tips, Curved, 4", Blades (104 mm)
- C** 744-250 STEVENS Tenotomy Scissors, 4-1/2" (11.4 cm), straight, sharp points
- 744-255 STEVENS Tenotomy Scissors, 4-1/2" (11.4 cm), straight, blunt points
- D** 744-268 STEVENS Tenotomy Scissors, 4-1/2" (11.4 cm), curved, sharp points
- 744-270 STEVENS Tenotomy Scissors, 4-1/2" (11.4 cm), curved, blunt point
- E** 744-280 AEBLI Corneal Scissors, 4" (10.2 cm), straight
- 744-282 AEBLI Corneal Scissors, 4" (10.2 cm), angled on flat, right
- 744-284 AEBLI Corneal Scissors, 4" (10.2 cm), angled on flat, left





A




B



C



D




E



F

RIBBON TYPE SCISSORS


- A 744-060 Eye Scissors, Ribbon Type, Pointed Tips, Straight, 4.1" Blade (22 mm)
- 744-062 Eye Scissors, Ribbon Type, Pointed Tips, Curved, 4.1", Blade (22 mm)
- B 744-064 STEVENS Tenotomy Scissors, Ribbon Type, Sharp Tips, Straight, 3.7", Blade (12mm)
- 744-066 STEVENS Tenotomy Scissors, Ribbon Type, Sharp Tips, Curved, 3.7", Blades (12 mm)
- C 744-068 STEVENS Tenotomy Scissors, Ribbon Type, Blunt Tips, Straight, 3.7", Blades (12 mm)
- 744-070 STEVENS Tenotomy Scissors, Ribbon Type, Blunt Tips, Curved, 3.7", Blades (12 mm)
- D 744-072 Strabismus Scissors, Ribbon Type, Straight, 4.2", Blades (25 mm)
- 744-074 Strabismus Scissors, Ribbon Type, Curved, 4.2", Blades (25 mm)
- E 744-076 Stitch Scissors, Ribbon Type, Extra Fine Blades, Light Curve, 3.7", Blades (20 mm)
- F 744-078 Utility Scissors, Ribbon Type, Blunt Tips, Straight, 4.1", Blades (20 mm)
- 744-080 Utility Scissors, Ribbon Type, Blunt Tips, Curved, 4.1", Blades (20 mm)



744-003




744-052



744-244

744-248

744-246



744-300

744-318



744-508


EYE INSTRUMENTS

- 744-003 Iris Scissors, Curved Sideways, 4-1/2" (11 cm)
- 744-052 Iris Scissors, Angular with probe point, 4-1/2" (11 cm)
- 744-244 BROPHY (Sullivan) Scissors, Straight, 5-1/2" (14.5 cm)
Sharp Points
- 744-246 BROPHY (Sullivan) Scissors, Curved, 5-1/2" (14.5 cm)
Sharp Points
- 744-248 BROPHY (Sullivan) Scissors, Curved, 5-1/2" (14.5 cm)
One Serrated Blade.
- 744-300 Strabismus Scissors, Straight, 4-1/2" (11 cm)
- 744-318 Strabismus Scissors, Curved, 4-1/2" (11 cm)
- 744-320 WESTCOTT Utility Scissors, Curved, Blunt Tips, 4-1/2" (11 cm)
- 744-508 BOWMAN Scissors, Curved, Probe Points, 4-1/2" (11 cm)




744-609

744-607



Castroviejo

744-600




Castroviejo

744-602




Vannas

744-994



Vannas

744-995



Vannas

744-996

744-566
744-567

744-557




744-565

A

EYE INSTRUMENTS

- 744-557 Enucleation Scissors, Straight, 5" (12.7cm)
- 744-565 Enucleation Scissors, light curve, 5" (12.7cm)
- 744-566 Enucleation Scissors, Medium curve, 5" (12.7cm)
- 744-567 Enucleation Scissors, full curve, 5" (12.7cm)
- 744-600 CASTROVIEJO Scissors, Curved, Sharp Points, 4" (10cm)
- 744-602 CASTROVIEJO Scissors, Curved, Blunt Points, 4" (10cm)
- 744-607 CASTROVIEJO Corneal Scissors, 3-3/4" (9.5 cm), angled blades 11 mm, blunt points
- 744-609 CASTROVIEJO Keratoplasty Scissors, 3-3/4" (9.5 cm), angled blades 11 mm, sharp points
- 744-994 VANNAS Micro Scissors, 3-1/4" (8.3 cm), straight
- 744-995 VANNAS Micro Scissors, 3-1/4" (8.3 cm), Curved
- 744-996 VANNAS Micro Scissors, 3-1/4" (8.3 cm), Angled on Flat



A

744-706



744-708



744-710



744-712



B

103-428



103-430



C

744-990



D

744-991



E

744-992

OPERATING SCISSORS

- 744-706 NOYES Iris Scissors, 4-1/2" (11.4 cm), straight, sharp/sharp points
- 744-708 NOYES Iris Scissors, 4-1/2" (11.4 cm), straight, sharp/blunt points
- 744-710 NOYES Iris Scissors, 4-1/2", straight (11.4 cm), blunt/blunt points
- 744-712 NOYES Iris Scissors, 4-1/2" (11.4 cm), curved, sharp/sharp points

- 103-428 BARRAQUER Iris Scissors, 2-1/4" (5.7 cm), angled 7 mm blades, sharp points
- 103-430 BARRAQUER Iris Scissors, 2-1/4" (5.7 cm), angled 7 mm blades, blunt points

- 744-990 MCPHERSON-VANNAS Micro Iris Scissors, 3-1/4" (8.3 cm), straight, sharp points
- 744-991 MCPHERSON-VANNAS Micro Iris Scissors, 3-1/4" (8.3 cm), Curved, sharp points
- 744-992 MCPHERSON-VANNAS Micro Iris Scissors, 3-1/4" (8.3 cm), Angled on Flat, sharp points




EYE INSTRUMENTS

- 745-000 VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), straight, extra delicate sharp points
- 745-018 VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), curved, extra delicate sharp points
- 745-020 VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), angled to side, extra delicate sharp points
- 745-026 VANNAS Capsulotomy Scissors, 3-1/4" (8.3 cm), angled on flat, extra delicate sharp points

- 745-028 Micro Scissors, 3-1/4" (8.3 cm), straight super fine blades
- 745-030 Micro Scissors, 3-1/4" (8.3 cm), curved super fine blades


- 745-042 MCGUIRE Corneal Scissors, 4-1/8" (10.5 cm), right
- 745-044 MCGUIRE Corneal Scissors, 4-1/8" (10.5 cm), left





OPERATING SCISSORS

- | | | |
|---|---------|---|
| A | 745-032 | Noyes Iridectomy Scissors 5 ½" (14 cm), Straight 7mm Blades |
| B | 745-035 | Castroviejo Corneal Section Scissors 4 ½" (11.4 cm), Right Inner Blade 1 mm Longer, Blunt Tips |
| | 745-040 | Castroviejo Corneal Section Scissors 4 ½" (11.4 cm), Left Inner Blade 1 mm Longer, Blunt Tips |
| C | 745-045 | McClure Iris Scissors 3 ¾" (9.5 cm), Angled Over Edge, 4mm Blades with Very Delicate Sharp Points |
| | 745-050 | McClure Iris Scissors 3 ¾" (9.5 cm) Angled on Flat, 4mm Blades with Very Delicate Sharp Points |
| | 745-055 | McClure Iris Scissors 4" (10.2 cm), Angled on Flat, Regular 10mm Blades with Sharp Points |



EYE INSTRUMENTS

- A** 746-562 Barraquer Needle Holder, 5", Curved w/out Lock
10mm Wide Hollow Round Handle
- 746-563 Barraquer Needle Holder, 5", Curved w/Lock
10mm Wide Hollow Round Handle
- 746-564 Barraquer Needle Holder, 5 1/8", Curved, Delicate Tapered Jaws
8mm Wide Round Handle
- B** 746-565 Kalt Needle Holder, 5-1/4", standard model, delicate jaws
- C** 746-566 Castroviejo Needle Holder, 5 1/2", Straight w/Catch
Extra Delicate Jaws
- 746-567 Castroviejo Needle Holder, 5 1/2" Curved w/Catch
Extra Delicate Jaws
- D** 320-432 Castroviejo Needle Holders, 5 1/4", Pen Type



EYE INSTRUMENTS

- A 749-507 Schiottz Tonometer, in Case, Standard
- 749-515 Schiottz Tonometer, in Case, Improved, Inclined Scale
- 749-560 Schiottz Tonometer, in Case, with Allen Modification, Inclined Scale

- B 749-564 Eye Magnet, with Loop and Magnifier, in Case
- C 749-606 Foreign Body Magnet with Probe and Loop

- D 749-903 Fox Eye Shields, Package of 12

- E 749-952 Castroviejo Caliper, 3-1/2"

- F 021-210 Flexible Stainless Ruler, 6" x 1/2" Graduated in 1/32" Fraction